

McLaren Print System Order

Order No: 7460
 Order Date: 2014-12-03
 User: Mike Stanton
 Phone: 810-342-5103

Ship Location: McLaren Flint - CCU 6th Floor
 401 S. Ballenger Hwy.
 Flint, MI

Brochures
 Quantity: 30
 Paragon Dept No: 20210
 Dept Name: McLaren Flint - CCU 6th Floor
 Company Number: 60

Order Total Price: 19.65

Item Number: M-24007
 Item Description: IV Insulin Adjustment Guide for Post Open Heart Patients Only
 Revision Date: 11/2014
 Print:
 Paper:
 Size:
 Fold:
 Finish:
 Drill:
 Misc Info: Size: 8.5 x 11 with lamination; 100 lb matte cover; CLC; SS

McLaren Flint IV Insulin Adjustment Guide

For Post Open Heart Patients Only
 Blood Glucose Goal 110-150 mg/dL

Insulin Titration Adjustment		CURRENT BLOOD GLUCOSE LEVEL (mg/dL)										
		less than 50	50-69	70-99	110-129	130-150	151-199	200-249	250-299	300-349	350-400	Greater than 400
PREVIOUS BLOOD GLUCOSE LEVEL (mg/dL)	less than 100			Hold	When insulin on HCLD. Restart insulin drip when BG greater than 110 mg/dL: - at 80% of previous rate. If held less than 4 hours - if held more than 4 hours, follow initiation table of signed orders							Give 10 units of NPH
	110-129			Hold	Continue	Continue	Increase by 20%	Increase by 20%	Increase by 50%	Increase by 50%	Increase by 50%	AND
	130-150	Hold & Give 10 units BID for 15-30 minutes until Blood Glucose less than 70 mg/dL	Hold & Give 10 units BID for 15-30 minutes until Blood Glucose less than 70 mg/dL	Hold	Decrease by 20%	Continue	Increase by 20%	Increase by 20%	Increase by 50%	Increase by 50%	Increase by 50%	
	151-199			Hold	Decrease by 50%	Continue	Increase by 20%	Increase by 20%	Increase by 50%	Increase by 50%	Increase by 50%	Increase rate by 50%
	200-249			Hold	Decrease by 50%	Decrease by 20%	Continue	Increase by 20%	Increase by 50%	Increase by 50%	Increase by 50%	
	250-299			Hold	Hold	Decrease by 50%	Decrease by 20%	Continue	Increase by 20%	Increase by 50%	Increase by 50%	Call Physician if Blood Glucose greater than 400 x 2
	301-349			Hold	Hold	Decrease by 50%	Decrease by 50%	Continue	Increase by 20%	Increase by 50%	Increase by 50%	
	350-400			Hold	Hold	Decrease by 50%	Decrease by 50%	Decrease by 20%	Continue	Increase by 20%	Increase by 50%	Increase by 50%
greater than 400			Hold	Hold	Hold	Decrease by 50%	Decrease by 50%	Decrease by 20%	Continue	Increase by 50%	Increase by 50%	

MAXIMUM INCREASE IN ONE HOUR IS 4 UNITS per HOUR / NO LIMIT ON RATE REDUCTION.

If the insulin drip is on HCLD, continue to check BG hourly.
 If on hold less than 4 hours, re-start the insulin infusion at 50% of previous rate when BG is greater than 110 mg/dL.
 If on hold more than 4 hours and BG greater than 110 mg/dL, Re-start insulin per insulin initiation chart on original signed insulin orders.
 Not to be used for DKA